



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Ceiling Shadowless OT Light with Dia 405mm X 2 Domes

**Product Code :**  
OT Equipment0035



**Description :**

Ceiling Shadowless OT Light with Dia 405mm X 2 Domes Key Structure The light is equipped with the reserve light source, which is activated in case the main light source brakes down. Lights arms are made of Stainless Steel Fitting. Sterilizable handle is provided in the middle of the dome. Inbuilt Parabolic Reflector, which increases the Lux up to 1,20,000 (10%). Heat absorbing & color correcting filter glass. Technical specifications Dome head : 405mm x 2 Halogen Bulbs : 24V x 150W x 4 nos. Light Intensity at 1 mt : 1,20,000 Lux Height Adjustment : 1014mm Parabolic reflector : 355mm Color Temperature : 4200 K 300 Temp rise in field : 2 - 3 C Control Unit : Transformer



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,