



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Infrared Ear Thermometer

**Product Code :**  
NMDIG10023



**Description :**

Infrared Ear Thermometer: Display range: 34C~44.0C(93.2F~111.2F), Accuracy: 0.2C(0.4F), Min Scale: 0.1, Memory: 10 memories, Measuring time: One Second, Beep function, Auto shut-off in 1 min, Battery: AAAx2, Alkaline batteries, Battery life: Approx 4000 times, Safe and hygienic Infrared Ear Thermometer, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers



**Medical Hospital Equipments Manufacturers**

**Website:** Medical Hospital Equipments Manufacturers,