



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**

Wall Mounted Shadowless OT Light with Single Dome Dia 465mm, Operated by Remote Control

**Product Code :**

OT Equipment0043



**Description :**

Wall Mounted Shadowless OT Light with Single Dome Dia 465mm, Operated by Remote Control Key Structure A New superior quality O.T. Light executed by a Remote Control. This Light provides the superior performance and is good for all kinds of Surgery. Lights arms are made of Stainless Steel Fitting. The light is equipped with the reserve light source, which is activated automatically in case the main light source brakes down. Sterilizable handle has been provided in the middle of the dome. LED Display that helps the user to adjust the Intensity of Light from 0 to 100% mode system. Technical specifications Dome head : 465mm x 1 Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 90,000 Lux Height Adjustment : 1078mm Dichoric Glass Reflector : 390mm Color Temperature : 4200 K 300 Temp rise in field : 2 - 3 C Control Unit : SMPS



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,