

**Product Name :**  
Baby Nipple Thermometer**Product Code :**  
NMBCE10014**Description :**

Baby Nipple Thermometer - Display range: 34C~44.0C(93.2F~111.2F, Accuracy: .1C(0.2F), Min Scale: 0.1, Fast, accurate and easy operation, Beep function, Digital display, Baby safe, Mercury free, Battery: 1.5V button battery(LR/SR<sup>-</sup>41), Battery life: about 100 hours under continuous operation. Baby Nipple Thermometer, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,