



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Baby Nipple Thermometer

Product Code :
NMBCE10014



Description :

Baby Nipple Thermometer - Display range: 34C~44.0C(93.2F~111.2F, Accuracy: .1C(0.2F), Min Scale: 0.1, Fast, accurate and easy operation, Beep function, Digital display, Baby safe, Mercury free, Battery: 1.5V button battery(LR/SR⁻41), Battery life: about 100 hours under continuous operation. Baby Nipple Thermometer, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers



Medical Hospital Equipments Manufacturers

Website: [SITE](#)