



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Trial Frame

Product Code :
OPHTH2010002



Description :

Hospital Trial Frame - 20 MM Reduced aperture. Hold 4 + 48 Lenses, I.P.D. Adjustment 50 mm to 78 mm Trial Frame, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers

Technical Specification :

Perfect for basic refraction needs, testing children or a practice back-up.
PD adjustable using sliding bridge mechanism.
The most simple and lightweight trial frame.
Light Weight And Fully Compact.
Holds 8 trial lenses.



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,