



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Halogen Headlight

**Product Code :**  
OPHTH2010008



**Description :**

Halogen Headlight: Light for true tissue color and consistent. Long lasting illumination. Iris diaphragm regulates diameter of light spot. The halogen illumination of the lamp offer an optimal light intensity which may be regulated by a switch. Halogen Headlight, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: SITE,