



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Halogen Headlight

Product Code :
OPHTH2010008



Description :

Halogen Headlight: Light for true tissue color and consistent. Long lasting illumination. Iris diaphragm regulates diameter of light spot. The halogen illumination of the lamp offer an optimal light intensity which may be regulated by a switch. Halogen Headlight, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,