



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Aspherical Lens

Product Code :
OPHTH2010011



Description :

14D provides high magnification along with an optimized working distance from the cornea, 20D diagnostic/laser lenses is perfectly corrected for field curvature, astigmatic aberrations, 30D is designed to provide the practitioner a wide field of view with a short working distance. Aspherical Lens, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,