



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Motorized Table

Product Code :
OPHTH2010017



Description :

Motorized Table: With switch mounted. Up & Down movement of the motorized table is aided by feather touch switch mounted with wooden top. This table has big lifting range with minimum height of 675 mm and a maximum of 875 mm. Power Supply: 220V or 110V, 50/60HZ, Up & Down movement by feather touch switch mounted with wooden top. Motorized Table, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,