



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Sintered Glassware

**Product Code :**  
NMCBS102174



**Description :**

Sintered Glassware We are indulged in offering premium quality Sintered Glassware which can be relied on for durability. Laboratory needs can be met with our impeccable range of Sintered Glassware. There are offered in variety of sizes and measurements. They are used for variety of purposes in lab. We are Manufacturer sintered glassware, gooch crucible, sintered disc, gas washing bottle supplier, salt bridge manufacturer, laboratory pipeline filter, conical funnel, buchner funnel with cone & stem, conical filter, buchner funnel, membrane filter holder assembly



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,