

**Product Name :**  
Wheelchair Folding (Travel)**Product Code :**  
NMWC0005**Description :**

Wheel chair folding, fix armrest, swing away footrest (pushing type) can easily be stored & can be carried in small car. 12.5cm castors in front & 31cm wheel at rear with brakes. This economical model is specially designed for traveling purpose. Wheelchair Folding (Travel), Wheel Chairs Exporters, Wheel Chairs Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,