

Product Name :
Saline Stands**Product Code :**
NMCBS90001**Description :**

Saline Stand To fulfill the needs and demands of our clients, we are able to provide Saline Stand that are available in different sizes and designs. All our products are widely appreciated by the customers for their corrosion & abrasion resistance, dimensional accuracy, robust construction and low maintenance. These products can also be customized as per the demands and needs of customers. Other Information: $1\frac{1}{2}$ HT. adjustable from Min.: 1370 mm to Max.: 2435 mm $1\frac{1}{2}$ Available in full stainless steel / epoxy powder coated stand & Stainless steel saline rod Saline Stand, Hospital Instrument exporters, Hospital Instrument manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:** Medical Hospital Equipments Manufacturers,