

**Product Name :**  
Hospital Bed With Built-In Backrest**Product Code :**  
NMHOPB00002**Description :**

Hospital Bed With Built-In Backrest from China Hospital Medical Furniture Manufacturers

**Technical Specification :**

**Hospital Bed With Built-In Backrest** Size : 198 L x 90 W cm x 60 H cms - Frame work made of strong rectangular steep tubes. Two section of CRCA M.S. Sheet uniformly perforated. Backrest section adjustable on ratchet mechanism. Head and leg bows of steel tubular pipes assembled to laminated pannels. Legs of ERW tubular peipes mounted on four heavy duty rubber shoes. Pre-treated and epoxy powder coated.

**Optional Accessories:** Side Railing, Saline rod, Mosquito Net Poles, Bed Mattress

Note : Above Bed can also supplied with castors.