

Product Name :
ORBITAL SHAKING INCUBATOR**Product Code :**
NMCBS11007**Description :**

ORBITAL SHAKING INCUBATOR Standard double wall construction. Outer chamber made out of mild steel finished in white stoving enamel/powder coated paint with mat finished colour combinations. Inner chamber made of stainless steel sheet. The double door is provided with double glass full view window. Interchangeable platform has been provided. OPTIONAL : i) Timer 0-24 hours. ii) Imported Thermostat iii) digital temp. Indicator-cum-Controller. Iv) Microprocessor PID digital temp. Indicator-cum-controller. Specification : ORBITAL SHAKING INCUBATOR Temperature Range : 5oC above ambient to 60oC + 1oC Speed Control : Upto 200 rpm Power Supply : 220V Ac 50Hz single phase Capacity : 16x500 ml conical flasks Make : Ajanta ORBITAL SHAKING INCUBATOR, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,