



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
WALK-IN-INCUBATOR

Product Code :
NMCBS11008



Description :

WALK-IN-INCUBATOR Designed for vast number of incubation applications. Can be installed even in hallway room or wherever convenient. Design feature includes enameled mild steel sheet exterior, aluminium / stainless steel, interior ventilators, internal lighting, exhaust fan etc. Suitable insulated 300 mm sq. Inspection window at the door permits inside view with -out opening the door. Complete with panel. Pre fabricated at factory assembled at site. Complete with cord and plug to work on 220 V single phase 50 Hz power supply. Temperature range from 5 $\text{ }^{\circ}\text{C}$ above ambient to 60 $\text{ }^{\circ}\text{C}$ \pm 1 $\text{ }^{\circ}\text{C}$. Controlled by digital indicator -cum -controller. Available sizes W x H x D a) 1.2x2.3x1.2 meter b) 2.3x2.3x2.3 meter WALK-IN-INCUBATOR, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,