

**Product Name :**  
SLIDE WARMING TABLE**Product Code :**  
NMCBS11016**Description :**

SLIDE WARMING TABLE Temperature range ambient temp. To 70° C i.e. Below the melting points of wax. The temperature is thermostatically controlled. The body is made out of thick mild steel duly finished in white stoving enamel/powder coated paint with mat finished colour combinations and top is made of highly polished aluminium sheet. To work on 220 V 50 Hz single phase Top plate size (Length x Width) a) 30x25 cm (12" x 10") b) 40x25 cm (16" x 10") c) 45x30 cm (18" x 12") d) 60x15 cm (24" x 6") SLIDE WARMING TABLE, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,