



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
SLIDE WARMING TABLE

**Product Code :**  
NMCBS11016



**Description :**

SLIDE WARMING TABLE Temperature range ambient temp. To 70° C i.e. Below the melting points of wax. The temperature is thermostatically controlled. The body is the made out of thick mild steel duly finished in white stoving enamel/powder coated paint with mat finished colour combinations and top is made of highly polished aluminium sheet. To work on 220 V 50 Hz single phase Top plate size (Length x Width) a) 30x25 cm (12â€•x10â€•) b) 40x25 cm (16â€•x10â€•) c) 45x30 cm (18â€•x12â€•) d) 60x15 cm (24â€•x6â€•) SLIDE WARMING TABLE, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)