

**Product Name :**  
AUTOCLAVE**Product Code :**  
NMCBS11019**Description :**

AUTOCLAVE Suitable for general laboratory use and sterilization. Stainless steel or aluminum vertical body. Domed cover can be sealed by joint-less gasket which is clamped to the vessel by wing nuts by hinged bolts which engage slotted lugs on rim of the cover. Provided with dial pressure gauge. Spring loaded safety valve has threaded cp for pressure adjustments. Supplied complete with plug & cord and is suitable to work on 2 20V AC 50Hz single phase. Specification :- AUTOCLAVE Size : a) 300x300 mm [Seamless] b) 300x375 mm [Jointed] c) 300x500 mm [Jointed] Pressure range : 15 psi to 18 psi Make : Ajanta AUTOCLAVE , Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers