



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
MORTUARY CHAMBER / FREEZER

Product Code :
NMCBS11022



Description :

MORTUARY CHAMBER / FREEZER Used to storing cadaverous under low temperature cool condition in order to prevent decomposition. These units are double walled with outer made of M.S. sheet duly painted and inner wall of stainless steel. Front opening covered door with magnetic gasket lock & key and handle individual chambers are provided independent lock with duplicate keys. Trays made of G.I. sheet to carry dead bodies are provided with handles. Trays travel on rollers fitted with ball bearing for easy movement. Hermetically sealed compressor serves refrigeration unit, to maintain temperature range of 2° C to 5° C or (-) 15° C to (-) 20° C (as ordered) Unit is fitted with Digital Temperature Indicator-cum-controller. Works on 220/230 volts A.C. Power supply. Capacity Temp. Temp. Voltage Stabilizer a) Single Body 2° to 5° C -15° to -20° C 3 KVA b) Two Bodies 2° to 5° C -15° to -20° C 4 KVA c) Three Bodies 2° to 5° C -15° to -20° C 5 KVA d) Four Bodies 2° to 5° C -15° to -20° C 8 KVA e) Six Bodies 2° to 5° C -15° to -20° C 8 KVA MORTUARY CHAMBER / FREEZER, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,