

**Product Name :**  
BLOOD BANK REFRIGERATOR**Product Code :**  
NMCBS11023**Description :**

BLOOD BANK REFRIGERATOR Blood Bank chamber is a most versatile unit for Blood Bank mounted on a sturdy angle frame. The inner chamber is made up of stainless steel sheet and outer of mild steel sheet duly powder coated paint. Caster wheels provided for easily mobility. The temperature is maintained in 4° C to 6° C with forced air circulation ensures uniform temperature, for storing blood bottles or blood bags. Temperature is controlled by Microprocessor PID digital temperature controlled-cum-indicator and triggers the AUDIO VISUAL ALARM in the event of fluctuation exceeding permitted tolerance. Unit is also fitted with DOOR OPERATED ALARM and POWER SUPPLY FAILURE ALARM, as a safety devices. A highly efficient refrigeration system consists of hermetically sealed compressor, for high performance cooling coil runs all around the cabinet to maintain uniformity of temperature. The entire refrigeration unit is housed at bottom. Two doors are provided. Inner door is of transparent acrylic sheet for easy viewing inside the chamber. Outer door is duly insulated for proper sealing from external environment and is lock able. Drawers provided for holding bags. Suitable BLOOD BANK REFRIGERATOR, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,