

Product Name :
Stability Chambers

Product Code :
NMCBS11028



Description :

Stability Chambers We offer a wide range of stability chambers that are precisely designed to create test chambers for carrying out studies on various environmental parameters such as humidity and temperature. These are also used in various application tests involving plants, tissues, microorganism and electronic components. These are available in two models STD (Outer M.S & Inner S.S) and GMP model (complete S.S.). Some of the specifications are: • Double walled construction inner SS – outer MS • 3 gap puff insulation • Full length inner glass door for clear inner vies • Outer metal door with magnetic gaskets & lock • Forced air circulation for uniform temperature throughout the chamber • Dual display preset electronic digital temperature controller cum indicator • Cooling system by means of hermetically sealed compressor unit coupled with evaporation coil and condenser • Temp. Range 10°C to 60°C 1°C • Humidity range 40% Rh to 98% Rh 3% Rh Our range is available in following models: Order Code Capacity Cu.ft Ltrs Inner Dimensions H x D x W (cms) No. of Shelves AEI-3a 3 cft 90Ltrs 50 x 43 x 43 2 AEI-3b 6 cft 180Ltrs 70 x 51 x 51 2 AEI-3c 10 cft 300Ltrs 100 x 55 x 55 3 AEI-3d 12 cft 360Ltrs 100 x 60 x 60 3 AEI-3e 16 cft 480Ltrs 114 x 65 x 65 4 Stability Chambers, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,