



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Cold Room

**Product Code :**  
MEDICAL-ANTM10010005



**Description :**

Cold Room

**Technical Specification :**

Cold room has total body with heavy sturdy angle iron & channels outer body is made of MS sheet duly painted. Inner chamber made S.S. single door opening system with door opening alarm. Temperature 4Â° C is brought down by hermetically sealed compressor and controlled by microprocessor based temperature variation alarm facility controller, puff insulation to maintain the temperature, inner chamber illumination total unit is moveable on casters.



**Medical Hospital Equipments Manufacturers**

Website: SITE,