

**Product Name :**  
Cadaverous Injector**Product Code :**  
MEDICAL-ANTM10010006**Description :**

Cadaverous Injector

**Technical Specification :**

Cadaverous Injector is used for injecting Formaldehyde solution in cadaverous at much higher speed than normal gravity process. Unit is fully covered & mounted on a portable trolley having four castor wheels for easy movement. Unit consists of one air compressor fitted with 1/2 HP Motor which is connected with a Stainless Steel Tank of 10liters capacity meant for storing and injecting the solution. Tank is fitted with a safety valve, pressure gauge and rubber tubing having provision for injection. Supplied completely with electric cord, plug and suitable to work on 220V AC 50Hz supply.