

Product Name :
Half Body CPR Emergency Skills Training Manikin

Product Code :
MED-EMG0004



Description :

Half Body CPR Emergency Skills Training Manikin

Technical Specification :

Simulating standard open airway
 External breast compression: display device and alarm device
 indicator light display of correct and wrong compression ; alarm of wrong compression;
 Intensity display of correct (at least 5cm) and wrong (less than 5cm) compression; alarm of wrong compression.
 Artificial respiration (inhalation): display device and alarm device
 Inhalation 600ml, wrong indicator light display and alarm prompting; inhalation between 500-600ml right indicator light display;
 Indicator light display open airway;
 Inhalation too quickly or too much result in air entering into stomach; indicator light display and alarm prompting.
 Ratio of compression and artificial respiration: 30:2 (one or two person).
 Operating cycle: one cycle includes five times of 30:2 ratio of compression and artificial respiration.
 Operation frequency: at least 100 times per minute
 Operation methods: exercise operation
 Examination of pupil response: mydriasis and myosis
 Examination of carotid response: pinch the pressure ball by hand and simulate carotid pulse
 Working conditions: Input power is 110-240V

Medical Hospital Equipments Manufacturers

Website: [SITE](#).

NAUGRA[®]