

**Product Name :**  
Half Body CPR Emergency Skills Training Manikin

**Product Code :**  
MED-EMG0004



**Description :**

Half Body CPR Emergency Skills Training Manikin

**Technical Specification :**

Simulating standard open airway  
 External breast compression: display device and alarm device  
 indicator light display of correct and wrong compression ; alarm of wrong compression;  
 Intensity display of correct (at least 5cm) and wrong (less than 5cm) compression; alarm of wrong compression.  
 Artificial respiration (inhalation): display device and alarm device  
 Inhalation 600ml, wrong indicator light display and alarm prompting; inhalation between 500-600ml right indicator light display;  
 Indicator light display open airway;  
 Inhalation too quickly or too much result in air entering into stomach; indicator light display and alarm prompting.  
 Ratio of compression and artificial respiration: 30:2 (one or two person).  
 Operating cycle: one cycle includes five times of 30:2 ratio of compression and artificial respiration.  
 Operation frequency: at least 100 times per minute  
 Operation methods: exercise operation  
 Examination of pupil response: mydriasis and myosis  
 Examination of carotid response: pinch the pressure ball by hand and simulate carotid pulse  
 Working conditions: Input power is 110-240V

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**