



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Industrial Refrigerators

**Product Code :**  
NMCBS11037



#### Description :

Industrial Refrigerators We offer a wide range of Refrigerators that provides an ideal freezing environment for advanced medical and industrial applications. These are basically used in prevention of plasma, related blood components, vaccines, micro-organism, and testing of materials and electronic components at -80 Deg C. Some of the features are: • Construction - Double Walled Construction, Exterior made of mild steel powder coated & interior made of stainless steel 304 (STD Model) or Exterior & interior made of stainless steel 304 (GMP Model)• High density CFC Free Poly Urethane Foam insulated cabinet• Keyed door lock for improved security with magnetic door gasket. • Full Length inner glass door for clear inner view. • Environment friendly Ultra Bio Freezers with micro control and CFC free cascade refrigeration system. • Air cooled, hermetically sealed compressor unit. • Temperature Range: -80°C controlled by microprocessor based temperature controllers. • Temperature alarm system with adjustable high - low alarm limits includes audio visual alarm on temperature deviation. • 230V Ac 50 HZ Mains. Further, the technical specifications are as follows: Order Code Capacity Cu.ft Ltrs Inner Dimensions H x D x W (cms) No. of ShelvesGS - 5a 3 cft 100L 50 x 45 x 45 2GS - 5b 6 cft 180L 70 x 51 x 51 2GS - 5c 10 cft 300L 100 x 55 x 55 3 Industrial Refrigerators, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)

