

Product Name :
Soxhlet Extraction Fat Analyzer

Product Code :
MED-PTH0008



Description :

Soxhlet Extraction Fat Analyzer

Technical Specification :

Small volume,fast heating speed,even heating.
Extraction bottle and apparatus are made of glass to avoid ether spills.
Accurate digital temperature controlling,easy to operate
Measuring Range: Oil content within the scope of 0.5%-60% of good,feed,oil and various kinds of fat product
Measuring Sample Quantity: 6 simultaneously
Heating time: Within 10 min.
Recovery System: Automatic
Recovery Rate: ?80%
Temperature Controlling Range: Room temperature 0-100?
Power: AC 220V/50Hz,1000W