

**Product Name :**  
Soxhlet Extraction Fat Analyzer

**Product Code :**  
MED-PTH0008



**Description :**

Soxhlet Extraction Fat Analyzer

**Technical Specification :**

Small volume,fast heating speed,even heating.  
Extraction bottle and apparatus are made of glass to avoid ether spills.  
Accurate digital temperature controlling,easy to operate  
Measuring Range: Oil content within the scope of 0.5%-60% of good,feed,oil and various kinds of fat product  
Measuring Sample Quantity: 6 simultaneously  
Heating time: Within 10 min.  
Recovery System: Automatic  
Recovery Rate: ?80%  
Temperature Controlling Range: Room temperature 0-100?  
Power: AC 220V/50Hz,1000W