

**Product Name :**  
Pathology Lab Full Auto Analyzer**Product Code :**  
MED-PTH0015**Description :**

Pathology Lab Full Auto Analyzer

**Technical Specification :**

Full Auto Analyzer offered by us has the ability to boost the number of samples that can be processed in the chemical testing laboratory. This type of analyzer has flexible sample tray and possesses 24-hour refrigeration of reagent tray. With Bi-directional LIS interface, the analyzer is very user-friendly and ensures simple operations

**Product Details :**

Application: Hospital

Number of Chambers: 1

Screen: LCD

Throughput: 100 tests per hour

Throughput with ISE: Up to 300 tests per hour (K, Na, Cl, Li)

Refrigeration Time: 24-hour refrigeration for reagent tray

Interface: Bi-directional LIS interface

**TYPICAL USERS**

Small laboratories (used as the main instrument); medium-sized laboratories (used as a back-up)

**FEATURES**

100 tests per hour, up to 300 tests per hour with ISE (K, Na, Cl, Li)

24-hour refrigeration for reagent tray

Flexible reagent/sample tray

Independent mixing stirrer

---

Robust and user-friendly operation software  
Best cost-efficiency  
Pre-dilution and post-dilution for sample  
Bi-directional LIS interface



**Medical Hospital Equipments Manufacturers**

**Website:** Medical Hospital Equipments Manufacturers,