

Product Name :
Pathology Lab Semi Auto Analyzer**Product Code :**
MED-PTH0016**Description :**

Pathology Lab Semi Auto Analyzer

Technical Specification :

The Semi-Auto Analyzer helps in measuring chemical indexes by examining blood and other body fluid. Further, it is combined with other clinical information, to assist in evaluating organs function, diagnosing diseases, recognizing disease gene and determining the norm for future therapy. The in-built thermal printer of the analyzer prints the results and makes, diagnosing a disease simpler process.

Product Details :

Application: Laboratory Use

Display Type: TFT touch-screen

Display Size: 7.0 inch

Number of Test Modes: 2

Test Mode Type: Flowcell or cuvette

Packaging Type: Boxes

SALIENT FEATURES

7.0" TFT touch-screen & pop-up keypad

Two test mode: flowcell or cuvette

Supports bi-chromatic tests for end point, fixed-time, kinetics methods

External keyboard or mouse support via USB

User-friendly interface and multi-language operation software

Built-in thermal printer



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,