


Product Name :
Pathology Lab Blood Cell Counter

Product Code :
MED-PTH0019



Description :

Pathology Lab Blood Cell Counter

Technical Specification :

Blood Cell Counter is completely automatic hematology analyzer with 19 parameters for microsampling technology and CBC test. Compact in size, the counter is generally used for counting the RBC in the blood. Further, the user-friendly interface of the counter allows simple and effective workflow. This type of counter ensures accurate results and efficient performance.

Product Details :

Application: Hospital, Clinical

Operation Mode: Semi Automatic

Sample Volume: 13 ul (whole blood)

Storage Up to: 20,000 sample results

Throughput: 30 samples per hour

Number of Counting Modes: 2

Counting Mode Type: Whole blood and Prediluted

FEATURES

Three-part differentiation of WBC, 19 parameters +3 histograms

Large Color LCD Display



Medical Hospital Equipments Manufacturers

Website: [SITE](#)