

**Product Name :**  
Pathology Lab Microscope

**Product Code :**  
MED-PTH0020



**Description :**

Pathology Lab Microscope

**Technical Specification :**

Microscope offered by us has quadruple revolving nosepiece and allows interchangeability of optics. Helps to clearly see the micro organisms, this is provided with mechanical stage right hand drive. The in-built 6v, 200W halogen light illuminator, in the microscope ensures high performance.

**Product Details :**

Application: Clinical, Laboratory

Surface Illumination: Halogen

Halogen Voltage Required: 6 V

Halogen Wattage: 200 W

Stage Length: 120 mm

Stage Width: 132 mm

**FEATURES**

Antifungus optics (sealed optical parts)

Interchangeability of optics

Abbe condenser & light relay system are fitted with high performance aspheric lenses for bright & crisp image  
Illumination system through SMPS circuit etc.

Coaxial coarse and fine roller guide (Comfortable position of focusing knobs). Tension adjustment on right

Quadruple revolving nosepiece

With mechanical stage right hand drive



## Medical Hospital Equipments Manufacturers

**Website:** Medical Hospital Equipments Manufacturers,