


Product Name :
Pathology Lab Centrifuge**Product Code :**
MED-PTH0021**Description :**

Pathology Lab Centrifuge

Technical Specification :

The Centrifuge operates under an irregular filtration process and is perfect to the competent separation of pharmaceutical products and bulk chemicals. This type of centrifuge is used along with tiny batch sizes and where products are changed often with less process costs. The digital speed indicator and Stepless speed regulator, makes this product highly efficient.

Product Details :

Application: Hospital, Clinical

Basket Material: SS

Display Type: Digital

Frequency: 50 / 60 Hz

Power Source: Battery

Packaging Type: Boxes

SALIENT FEATURES

Digital Pace Indicator

Stepless speed controller

0-60 minutes digital countdown timer

Press Switch for quick acceleration to full speed



Medical Hospital Equipments Manufacturers

Website: [SITE](#)