

Product Name :
Pathology Lab HB1Ac Machine**Product Code :**
MED-PTH0036**Description :**

Pathology Lab HB1Ac Machine

Technical Specification :

Sample type: Venous blood, Peripheral finger blood
Sample volume: 5?l
Sampling mode: Full-auto 25-position sample turntable
Testing mode: Auto test for batch
Testing item: Glycated Hemoglobin HbA1c
Display: 320x240 LCD Graphic Display
Data storage: 1000 test report
Power supply: AC220V±22V 50Hz 65VA