

Product Name :
Pathology Lab Eliza Plate Washer and Reader**Product Code :**
MED-PTH0037**Description :**

Pathology Lab Eliza Plate Washer and Reader

Technical Specification :

8-channel vertical optical system.

Super fast reading speed of 7 seconds.

10.1-inch Large Touch Screen with visual layout function, what you see is what you get.

Hardware: ARM architecture high-performance low-power processor, with ADI high-precision AD conversion chip, the measurement is more accurate and stable.

Software: Self-developed core algorithm, with independent intellectual property rights, flexible use of multiple calculation methods.