

**Product Name :**  
Pathology Lab Immunoassay Analyzer

**Product Code :**  
MED-PTH0038



**Description :**

Pathology Lab Immunoassay Analyzer

**Technical Specification :**

Modes of Operation: Random, Batch and STAT  
Separation method: Magnetic separation technology  
Sample Positions: 60 (Each position could be used as an emergency position)  
Reagent Positions: 30 pcs (2-8°C refrigeration)  
Incubation Positions: 50  
Time to 1st Result: 18 Minutes  
Sample Volume: 15-100ul  
Detection Precision: CV?5%  
Carry-over Rate: ?10-5  
Correlation Coefficient: r?0.99  
Calibration Stability: 28 Days  
Calibration Type: 6/7-point calibration  
Software system: Windows 7/8/10, 32 or 64 bit  
Bi-direction, support: HL7 protocol  
Intuitive User Interface, Intelligent Data Management  
Interface: TCP/IP Network interface  
Other Function: Liquid level detection  
Anti-collision function  
Sample barcode scanning  
Reagent radio frequency identification  
Intelligent alarm prompt

---

Power Supply: AC110/220V,50/60Hz

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,