

Product Name :
Pathology Lab Grossing Workstation

Product Code :
MED-PTH0044



Description :

Pathology Lab Grossing Workstation

Technical Specification :

External Size (W*D*H): 1800*750*1950mm
Working Area Size (W*D*H): 1695*430*900mm
Work Surface Height: 900mm
Consumption: 300W
Optional water heater: 1500W
Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz
Exhaust Airflow Volume: 500-600m³/h
Material Main Body: 304 stainless steel
Work Table: 304 stainless steel
Standard Accessories: Spotlight;
UV lamp;
Fluorescent lamp;
Exhaust system;
4m exhaust duct;
Pipe strap;
Auto flushing device;
Stainless steel hanger;
Adjustable hot and cold water tap;
180mm deep stainless steel water sink, funnel with filter.
Active carbon filter
Optional Accessories: Pulverizer, Electrical water heater, Automatic induction hand dryer, Eye washer, Camera.

Flex magnifier light, laptop bracket

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,