

**Product Name :**  
Pathology Lab Tissue Floatation Bath**Product Code :**  
MED-PTH0045**Description :**

Pathology Lab Tissue Floatation Bath

**Technical Specification :**

Temp. Range: RT to 90°C

Temp. Precision: 1°C

Environmental Temp. Range: 0~40°C

Bath Bowl Size(W\*D\*H): 240\*180\*50mm

Display: LED

Power Supply: AC220V±10%, 50/60Hz; AC110V±10%, 60Hz

Consumption: 200W