



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Water Baths Constant Temperature (Precision)

Product Code :
NMCBS11043



Description :

Water Baths Constant Temperature (Precision) We offer a wide range of water baths with constant temperature. Our range of water baths are including water bath (round), water baths round thermostatic, water baths (rectangular), water bath rectangular thermostatic, water baths universal and water bath serological etc. Some of the features are: •Rectangular in shape •Double walled –Inner SS – Outer MS/ Complete S.S •Temperature Range : 5°C above ambient to 95°C $\pm 0.1^\circ\text{C}$; •Water circulatory pump with stirrer-Precision model- •Digital Electronic Preset Temperature controller and Immersion Heaters. Water Baths Constant Temperature (Precision), Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,