

**Product Name :**  
Pathology Lab Tissue processor**Product Code :**  
MED-PTH0046**Description :**

Pathology Lab Tissue processor

**Technical Specification :**

Fully intelligent design, ensure timely automatic determination and recovery from an abnormal event.  
Manual adjustment can be conducted anytime during the programmed automatic operation; afterwards, the system automatically enters the programmed operation.

Green inner-cycling air purification system to efficiently absorb and remove poisonous gas.

Automated fan and light control.

Internal dry heating mechanism with high-precision temperature control

Cups Number: 12

Single Cup Capacity: 1L

Processing Time Within 0~99 hours for the first cup(BK-TS3D Any length); 0-24 hours for other cup.

Paraffin Cup Temp. Range: RT~99°C

Temp. Precision: 1°C

Dripping Time: 0~30 s 0~60s

Display: LCD

Agitation Frequency: 2 times/min

Environment Temp. Range: 0~40°C

Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**