

**Product Name :**  
Pathology Lab Water Bath**Product Code :**  
MED-PTH0048**Description :**

Pathology Lab Water Bath

**Technical Specification :**

Stainless steel inner chamber and upper cover.  
Adopt new design structure, which could avoid the leakage of steam.  
PID Microprocessor controller, LED display.  
Capacity: 3L  
Temp. Range: RT.-100?  
Temp. Fluctuation:  $\pm 0.5$ ?  
Temp. Accuracy:  $\pm 1$ ?  
Time Range: 0~999min  
Power: 500W  
Voltage: AC110/220VV $\pm 10\%$ , 50/60Hz