

**Product Name :**  
Pathology Lab Humidity Cabinet

**Product Code :**  
MED-PTH0052



**Description :**

Pathology Lab Humidity Cabinet

**Technical Specification :**

External Size(W\*D\*H): 700\*650\*1230 mm  
 Internal Size(W\*D\*H): 600\*500\*540 mm  
 Tested Opening: Safety Height ? 200 mm(8")  
 Max Opening: 370mm(15")  
 Inflow Velocity: 0.53±0.025 m/s  
 Down Flow Velocity: 0.33±0.025 m/s  
 HEPA Filter: Two, 99.999% efficient at 0.3µm. Filter life indicator.  
 Front Window: Motorized. Toughened glass. Anti UV.  
 Noise: EN12469 ? 58 dB / ? 61 dB  
 UV Lamp 15W\*1: UV timer, UV life indicator, emission of 253.7 nanometers for most efficient decontamination.  
 LED Lamp: 8W\*2  
 Illumination: ?1000Lux  
 Consumption: 500W  
 Socket: Two, total load of two sockets: 500W  
 Display: LCD display  
 Control System: Microprocessor  
 Airflow System: 70% air recirculation, 30% air exhaust  
 Visual and Audio alarm: Filter replacement, window over height, abnormal air flow velocity.  
 Material: Work Zone: 304 stainless steel  
 Main Body: Cold-rolled steel with anti-bacteria powder coating.

---

Work Surface Height: 750mm with optional base stand  
Power Supply: AC220V $\pm$ 10%, 50/60Hz; 110V $\pm$ 10%, 60Hz

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,