

**Product Name :**  
Pathology Lab Freezer Refrigerator

**Product Code :**  
MED-PTH0055



**Description :**

Pathology Lab Freezer Refrigerator

**Technical Specification :**

External Size(W\*D\*H) mm: 500\*550\*1300  
Capacity: 130L  
Temp. Range: 2°C~8°C  
Temp. Accuracy: 0.1°C  
Control System: Microprocessor Control, Large Screen LED Display  
Alarm: Audible and visual alarm for: High and low temperature, Power failure alarm, Sensor failure, Door ajar  
Refrigeration Type: Forced air refrigeration system  
Refrigerant: R600a, CFC Free  
Condenser & Evaporator: Bundy tube condenser  
Defrost: Auto defrost  
Door: Glass door with heat reflection film  
Construction: Structure Unibody Design & Mono-assembly foaming (Rigid polyurethane insulation material)  
Internal  
Material: Embossed aluminum sheet  
External  
Material: PCM  
Shelves 4 pcs  
Door: Glass door with heat reflection film  
Consumption: 136W  
Power Supply: AC110/220V±10%, 50/60Hz



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,