

Product Name :
Pathology Lab Mixers**Product Code :**
MED-PTH0058**Description :**

Pathology Lab Mixers

Technical Specification :

Shaking Mode: Orbital
Orbital Diameter: 4mm
Motor Input: 58W
Motor Output: 10W
Speed Range: 0~2500rpm
Speed Display: Scale
Operation Mode: Continuous/touch operation
Ambient Temp: 5-40°C
Relative Humidity: 80%
Protection Type: IP21
Power Consumption: 60W
Power supply: AC110/220V±10%, 50/60Hz

