

Product Name :
Pathology Lab Stirrers**Product Code :**
MED-PTH0060**Description :**

Pathology Lab Stirrers

Technical Specification :

Work Plate Size: 120*120mm

Work Plate Material: Aluminum

Max. Stirring Capacity: 2L

Speed Range: 100-2000rpm

Heating Temp: 380 °C

Power Consumption: 180W

Power Supply: AC110/220V±10%, 50/60Hz

Standard Accessory: Adjustable rack, stirrer bar; PT100 sensor (only for BS-4HC)