

**Product Name :**  
Serological Water Bath**Product Code :**  
NMCBS11047**Description :**

Serological Water Bath Available in two different models - exterior made of mild steel powder coated and interior made of stainless steel 304(std model) with piramidal ss lid & exterior made of stainless steel 304 and interior made of stainless steel 316 (gmp model) with piramidal ss lid. temperature range: 5deg. C above ambient to 85deg. C +/-2deg. C .

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,