



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ultra-Low Degree Deep Freezer Laboratory

Product Code :
NMCBS11058



Description :

Ultra-Low Degree Deep Freezer Laboratory - Vertical Deep Freezer for **Hospital Refrigeration System**
Exporters, **Medical Lab Refrigeration Equipments** Manufacturers

Technical Specification :

Ultra-Low Degree Deep Freezer Laboratory - **Vertical Deep Freezer Capacity 158L**
Condenser & Evaporator Made of copper

Construction Structure:

Unibody Design & Mono-assembly foaming (Rigid polyurethane insulation material)

Internal SS304

External Material Cold-rolled steel with anti-bacteria powder coating

Door: lockable door

Caster: 4 casters



Medical Hospital Equipments Manufacturers

Website: [SITE](#)

