



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Round & Rectangular Vacuum Oven S

**Product Code :**  
NMCBS11061



**Description :**

Vacuum Oven S (Round & Rectangular) We offer a wide range of vacuum ovens that are specially designed to provide excellent vacuum & temperature. These are manufactured from quality materials thus providing maximum satisfaction to our clients. We offer these ovens at most competitive prices. Some of the features are as follows: •Construction: Double walled construction, Interior made thick stainless steel 304 & Exterior made of mild steel powder coated finish. •Insulation: The gap between the walls is filled with high grade mineral glass wool. •Door: Aluminum casting door with silicone rubber gasket. Door fitted with Toughened glass window for inner observation. •Temperature control: Digital display temperature controller / Micro Processor based PID controller. •Temperature Range: 5°C above ambient to 130 °C/250 °C •Vacuum gauge & Cock provided. •0 to 760 mm of Hg vacuum. Vacuum Oven S (Round & Rectangular), Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,