



Email : sales@naugramedical.com

Phone: 9896600003

**Product Name :**  
Water Baths Constant Temperature (Precision)

**Product Code :**  
NMCBS11062



**Description :**

Water Baths Constant Temperature (Precision) We offer a wide range of water baths with constant temperature. Our range of water baths are including water bath (round), water baths round thermostatic , water baths (rectangular), water bath rectangular thermostatic , water baths universal and water bath serological etc. Some of the features are: •Rectangular in shape •Double walled –Inner SS – Outer MS/ Complete S.S •Temperature Range : 5°C above ambient to 95°C  $\pm 0.1^\circ\text{C}$ ; •Water circulatory pump with stirrer-Precision model- •Digital Electronic Preset Temperature controller and Immersion Heaters. Size Dimensions L x B x H (cms) Capacity a) 30 x 25 x 20 15 L b) 35 x 30 x 25 26 L c) 45 x 35 x 30 47 L Water Baths Constant Temperature (Precision), Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,